



10.36 Overbed table with castors



Categoria Prodotto: [Beds and Homecare Aids](#)

URL Prodotto: <https://vassilli.it/en/prodotto/10-36-overbed-table-with-castors/>

Descrizione Prodotto

Varnished steel frame that can be disassembled. Laminate top. Top with raised edges on the horizontal sides to avoid the fall down of objects once reclined.

Adjustment in height and inclination of the top. 4 twin wheels at the base two with brake.

Dim. cm 70x45x68/110h. Weight kg 10,5. Volume m3 0,05.

Requirements following Rule ISO EN 12182.