



## 10.58E Hydraulic base adjustable in height



**Categoria Prodotto:** [Beds and Homecare Aids](#)

**URL Prodotto:**

<https://vassilli.it/en/prodotto/10-58e-hydraulic-base-adjustable-in-height/>

### Descrizione Prodotto

Structure made of structural steel powder varnished. Frame with cross for the height adjustment from cm 40 to cm 80 by hydraulic pump by foot control on the 2 sides for up and down.

Connections for application to our beds.

**Possibility to insert: set of wheels Ø 125 art. 16.01.**

Dim. Cm 160x90x40+40h.

Load capacity kg 230.

Requirements following Rule CEI UNI EN 60601-2-52.

**Available also for beds with width cm 120 art. 10.58E/120.**