



10.72EN Electric bed adjustable in height



Categoria Prodotto: [Beds and Homecare Aids](#)

URL Prodotto:

<https://vassilli.it/en/prodotto/10-72en-electric-bed-adjustable-in-height/>

Descrizione Prodotto

The “Nursing” panels bordered in tubular steel, removable without tools, middle panels made of medium density laminate wooden finishing.

All the components used for the manufacture and finishing are atoxic and biocompatible. Possibility to insert: set of wheels art. 16.01, overbed trapeze bar applied to the bottom art. 10.44, I.V. pole art. 13.02L.

Upsetting and removable side rails art. 10.56-10.53.

Requirements following Rule CEI UNI EN 60601-2-52.

Dim. cm 214x90x40+45h net. Load capacity kg 230. Available also on size width cm 120 art. 10.72EN/120.