



10.72ENT Electric bed adjustable in height



Categoria Prodotto: [Beds and Homecare Aids](#)

URL Prodotto:

<https://vassilli.it/en/prodotto/10-72ent-electric-bed-adjustable-in-height/>

Descrizione Prodotto

Top panels, removable without tools setted up for insert the side rails art. 10.59 “top” and art. 10.59A aluminium side rails.

All the components used for the manufacture and finishing are atoxic and biocompatible.

Possibility to insert: set of wheels Ø 125 art. 16.01, overbed trapeze bar applied to the bottom art. 10.44, I.V. pole art. 13.02L.

Dim. Cm 217x104x40+45h net. Load capacity kg 230.

Requirements following Rule CEI UNI EN 60601-2-52.