



10.72TR Therapeutic bed



Categoria Prodotto: [Beds and Homecare Aids](#)

URL Prodotto: <https://vassilli.it/en/prodotto/10-72tr-therapeutic-bed/>

Descrizione Prodotto

New Therapeutic bed destined to home care for seriously ill patients and/or terminal

It could be an opportunity for the Health Systems and/or homecare assistance service the availability of a very technical aid having therapeutic value which improves the life of patients having very serious illness because in addition to avoid further pathologies and complications to the patient, it could be a way to reduce the costs connected to the unavoidable complications which would arise without using this aid.

Articulated 4 sections bottom with movements controlled by electric actuators. Base with 4 pivoting wheels Ø 125 mm with pedal brake. Head, foot panels and

side rails made of solid beech wood with walnut panels painted in natural. The backrest net during the inclination movement performs simultaneously a movement of horizontal rear transfer with articulation reciprocal to the deformation of the mattress, in order to maintain the user in a stable position avoiding the classic push pressure towards the footboard. Legs net reclining until 45° and foot area to follow by side pantographs. Bed bottom adjustable in height with electric movement operated by 2 telescopic independent columns placed at the edges of the bed.