



10.77CS200 COMPACT STAND-UP HOIST



Categoria Prodotto: [Extra Large Solutions](#)

URL Prodotto: <https://vassilli.it/en/prodotto/10-77cs200-compact-stand-up-hoist/>

Descrizione Prodotto

The **Compact Stand-Up hoist** has been designed to meet the needs of practicalness, safety, flexibility and comfort which are requested by the daily use of such products. The small size allows the passage in limited spaces and the full control by the attendants.

The components as the lifting arm adjustable in depth, the column adjustable in height, the footplate and the knee restrain and the slings adjustable, allow to make the lifting operations in full safety and harmony with the anthropometric features of the user.

Equipped with **corset slings with adjustment rings** available in sizes for different bodies. The legs of the base, equipped with sturdy twin wheels (2 with foot brake), can be opened by lever for easy frontal access to wheelchairs or

other in full stability and safety.

The movement of the arm is guaranteed by the electric actuator with telecommand with extensible cable. Control box with emergency switch and removable batteries. Equipped with battery charger. **Max load capacity kg. 200**