



## 10.77CSN200 COMPACT ACTIVITY LIFTER - Kg 200



Categoria Prodotto: [Transfer](#)

URL Prodotto: <https://vassilli.it/en/prodotto/10-77csn200-compact-activity-lifter-kg-200/>

### Descrizione Prodotto

**Compact activity lifter solution with load capacity 200 Kg.**

Technical and mechanical features as per **art. 10.77C**; additional features: the column adjustable in height, the adjustable footplates, the adjustable knee restrain and adjustable corset slings , allow to make the verticalizing in full safety and harmony with the anthropometric features of the user.

The reduced dimensions allow the passage on narrow passages and the full control by the assistants.

Equipped with corset slings with adjusting rings; in alternative are available different sizes for different builds. Equipped with automatic charger.

**15.99B**

Additional batteries with charger separate from the lifter, allows continuity of use.