



10.77N200 ELECTRIC BATTERY LIFTER



Categoria Prodotto: [Transfer](#)

URL Prodotto: <https://vassilli.it/en/prodotto/10-77n200-electric-battery-lifter/>

Descrizione Prodotto

The **patient lifter** is tested for a load until **150 kg**.

Its lifting stroke allows to perform in an easier way all the transfers which are necessary in the daily care. The adjustment to the patient mobility is obtained by various types of slings.

Thanks to the **twin castors 2 of them with brake**, it is possible to move easily under the beds usually used at home and at hospitals. The opening base allows **easy transfer** from and to the wheelchair.

Handset with 2 mt of cable, power module with emergency button and led for battery charge status, 2 batteries 12V. Possibility to insert all the models of slings of our programme.

Maximum load capacity 200 kg.