



## 10.77T TOWER LIFT



**Categoria Prodotto:** [Extra Large Solutions](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/10-77t-tower-lift/>

### Descrizione Prodotto

#### HOIST WITH ARTICULATED HANGER WITH LOAD CAPACITY UNTIL KG 280

The **Tower Lift** is the dynamic and complete solution for the moving of patients in safety, thanks to its stability. The column allows a stroke in height of 110 cm approximately. It is possible to pick up a patient from the floor. The hanger where are hooked the slings allows to rotate the patient but it doesn't tilt during the moving.

**Tower Lift** is a lifter with vertical electric moving, also the leg opening is electric, all is feeded on low tension at 24V by a telecommand complete with battery box, electronic system and battery charger and it satisfies all the safety standards forecasted by the Rule 10535 UNI EN ISO.

Slings on sizes: L-XL-XXL.

**Max load kg 280**

Available slings: ART. 15.94XXL Stretcher slings.