



10.77T200 TOWER LIFT



Categoria Prodotto: [Transfer](#)

URL Prodotto: <https://vassilli.it/en/prodotto/10-77t200-tower-lift/>

Descrizione Prodotto

The **Tower Lift** is the dynamic and complete solution for the moving of patients in safety, thanks to its stability. The column allows a stroke in height of 110 cm approximately. It is possible to pick up a patient from the floor. The hanger where are hooked the slings allows to rotate the patient but it doesn't tilt during the moving.

Tower Lift is a lifter with vertical electric moving, also the leg opening is electric, all is feeded on low tension at 24V by a telecommand complete with battery box, electronic system and battery charger and it satisfies all the safety standards forecasted by the Rule 10535 UNI EN ISO.

Slings on sizes: L-XL-XXL. Max load kg 200.

Available slings:

- ART. 15.95A Slings for bilateral amputated. - ART. 15.95A-RT/LF Slings for unilateral amputated. - ART. 15.95 Hammock slings with head support. - ART. 15.95E Slings with under axilla catch and slings for the legs. - ART. 15.94 Stretcher slings. - ART. 13.07T Dynamometer with hanger. - ART. 15.93T Articulated hanger for tower lift. - ART 15.95T Wrapping up sling for articulated hanger