



## 10.79M MOLEKOLA



Categoria Prodotto: [Transfer](#)

URL Prodotto: <https://vassilli.it/en/prodotto/10-79m-molekola/>

### Descrizione Prodotto

**Molekola** is a basic aid to help the **transfers of cooperative users** (not seriously ill) from the bed to the wheelchair and viceversa, from the wheelchair to the WC and viceversa, from the wheelchair to the car and viceversa, ecc.

The innovative aid has been designed by Tecnothon and made by Vassilli S.r.l. Italy, its total weight is less than 2 Kg. and can be used for persons having weight until 100 Kg. due to the high resistant materials used. It is reducible and can be inserted also in the suitcase to be used on travel or at holidays houses.

The dimensions once closed are 50x14cm, the central column is adjustable in height and can be used in different living situations. The finishing can be varnished in brilliant colours or in carbon or natural polished.