



## 11.27A Sliding Walker on 4 castors complete with antibrachial supports



**Categoria Prodotto:** [IHS](#)

**URL Prodotto:**

<https://vassilli.it/en/prodotto/11-27a-sliding-walker-on-4-castors-complete-with-antibrachial-supports/>

Descrizione Prodotto

Walker with 4 wheels Ø 100 mm at the base, the front ones steering. Soft padded axillary supports adjustable in height with lever clamp. Front soft rubber handgrips adjustable in depth, with hand lever brakes for emergency and parking acting on the rear wheels. Folding frame.

Load capacity kg 100.

Requirements following Rule EN ISO 11199-2.