



11.27T Sliding Walker on 4 castors complete with tray



Categoria Prodotto: [IHS](#)

URL Prodotto: <https://vassilli.it/en/prodotto/11-27t-sliding-walker-on-4-castors-complete-with-tray/>

Descrizione Prodotto

Padded chest laying tray covered in soft plastic material washable and antiscratch. Front soft rubber handgrips, adjustable in depth with hand lever brakes for emergency and parking acting on the rear wheels. Support structure of the tray adjustable in height with lever clamps. Folding frame.

Load capacity kg 100.

Requirements following Rule EN ISO 11199-2.