



11.40N Walker with antibrachial supports



Categoria Prodotto: [IHS](#)

URL Prodotto:

<https://vassilli.it/en/prodotto/11-40n-walker-with-antibrachial-supports/>

Descrizione Prodotto

Structure made of varnished steel, foldable, antibrachial PU supports adjustable,

wrapping, with handles and emergency+ parking adjustable brakes, polyurethane seat, basket, four puncture proof wheels Ø 200x50, the front ones driving type on ball bearings. Load capacity kg 100.

Requirements following Rule EN ISO 11199-2.