



11.40T Walker with padded tray



Categoria Prodotto: [IHS](#)

URL Prodotto: <https://vassilli.it/en/prodotto/11-40t-walker-with-padded-tray/>

Descrizione Prodotto

Structure made of varnished steel, foldable, chest padded laying tray covered in soft plastic material washable and antiscratch. Front soft rubber handgrips,

adjustable in height and depth, with hand lever brakes for emergency and parking which act on the rear wheels, polyurethane seat, basket, 4 puncture proof wheels Ø 200x50 mm. the front ones steering on ball bearings.

Load capacity kg 100.

Requirements following Rule EN ISO 11199-2.