



11.42N Walker type rollator mod. Fiesta



Categoria Prodotto: [IHS](#)

URL Prodotto:

<https://vassilli.it/en/prodotto/11-42n-walker-type-rollator-mod-fiesta/>

Descrizione Prodotto

Folding structure made of varnished aluminium, handles adjustable with soft

handgrips and emergency + parking brakes, padded seat, tubular padded backrest, basket with removable handle, wheels Ø 200x35 puncture proof, the two front ones revolving and driving on ball bearings.

Load capacity kg 100.

Requirements following Rule EN ISO 11199-2.