



11.50/2RG JUNIOR FIXED AND RECIPROCAL WALKER



Categoria Prodotto: [Junior Line](#)

URL Prodotto: <https://vassilli.it/en/prodotto/11-50-2rg-junior-fixed-and-reciprocal-walker/>

Descrizione Prodotto

Suitable for childrens from 8 to 15 years

Aluminium frame adjustable in height from cm 68 to cm 75. Possibility to articulate the frame disengaging the lockage which inhibits the articulation. Soft polyurethane handgrips. Two front revolving wheels and two rear tips. Folding frame. Weight kg 3,7. Load capacity kg 100.

Requirements following the Rule EN ISO 11199-2.