



11.50/2RG Walker



Categoria Prodotto: [IHS](#)

URL Prodotto: <https://vassilli.it/en/prodotto/11-50-2rg-walker/>

Descrizione Prodotto

Aluminium frame adjustable in height from cm 74 to cm 92. Possibility to articulate the frame disengaging the lockage which inhibits the articulation. Soft polyurethane handgrips. Two front revolving wheels and two rear tips. Folding frame. Weight kg 3,7.

Load capacity kg 100.

Requirements following Rule EN ISO 11199-2