



11.50/4RP 4 wheels sliding walker with pressure brake



Categoria Prodotto: [IHS](#)

URL Prodotto: <https://vassilli.it/en/prodotto/11-50-4rp-4-wheels-sliding-walker-with-pressure-brake/>

Descrizione Prodotto

Structure made of tubular polished duraluminium with anodizing protection treatment, with steel powder varnished components. Frame folding on the upper cross by quick release button. 4 telescopic legs for the height adjustment with pre holed locations for quick release. 2 front revolving wheels, 2 rear fix with handle pressure emergency and parking brake which by the effect of the load exercised insert the end cap at floor. Soft polyurethane shaped antislip handgrips. All the materials used for the manufacture of the aid are hypoallergenic and atoxic. Load capacity kg. 100.

Requirements following Rule EN ISO 11199-2