



11.50/4RP JUNIOR WALKER WITH 4 WHEELS WITH PRESSURE BRAKE



Categoria Prodotto: [Junior Line](#)

URL Prodotto: <https://vassilli.it/en/prodotto/11-50-4rp-junior-walker-with-4-wheels-with-pressure-brake/>

Descrizione Prodotto

Suitable for childrens from 8 to 15 years]

Aluminium frame folding and adjustable in height. Possibility to articulate the frame disengaging the lockage which inhibits the articulation. 2 front pivoting castors, 2 rear fix castors with emergency pressure brake on the upper handles which consequently to the load pressure exercised inserts the floor end tip. Soft polyurethane handgrips.

Requirements following the Rule EN ISO 11199-2. Load capacity kg 80.