



11.71AS ALL SERVICE



Categoria Prodotto: [Comode](#)

URL Prodotto: <https://vassilli.it/en/prodotto/11-71-as-all-service/>

Descrizione Prodotto

The aid is destined to moving on very narrow spaces like airplane corridors and/or trains and/or boats for passengers transport. The special structure of the frame allows the passage on corridors and/or strait spaces with a width of cm 40.

The wheelchair has been designed on request to overpass such situations. The wheelchair is fitted up with safety chest adjustable belt, emergency and parking brakes operating by the attendant, upsetting and lowering armrests to simplify the side transfer. Seat, backrest and armrests are padded and washable. The footplates are removable and made of composite material. Two front wheels Ø320 mm with brake, two rear wheels Ø180 mm all of puncture proof type.

By the particular structure and with the insertion of an oildynamic pump foot operating, the seating can be adjusted from 50 to 80 cm, in such way also a passenger of greater dimensions than the seating can have access on because the seating raises and can overcome eventual obstacles.