



## 11.71N Wheelchair for patienttransport



**Categoria Prodotto:** [Comode](#)

**URL Prodotto:**

<https://vassilli.it/en/prodotto/11-71n-wheelchair-for-patienttransport/>

### Descrizione Prodotto

Varnished steel frame. Synthetic, heavy duty, washable material upholsteries. Side adjustable and elevating footrests to allow the patient to seat down and raise up,

while avoiding the tipping of the wheelchair.

2 front wheels Ø 300 (11,81") with brakes, 2 rear ones Ø 200 (7,87"). Crush - proof polyurethane armrests upsetting with front lock.

Rear handgrips for the attendant. Backrest and armrest can be disassembled for space saving during transport. Fitted up with safety belt.