



11.86S COMODA



Categoria Prodotto: [Comode](#)

URL Prodotto: <https://vassilli.it/en/prodotto/11-86s-comoda/>

Descrizione Prodotto

Varnished steel frame. High density foam paddings. Upholstery without sewing, synthetic washable antiscratch fabric.

Shaped seat with W.C device complete with removable box. Removable side panels

with polyurethane armrests. Removable legrests with footplates opening and adjustable in height. 4 wheels Ø100mm, the rear ones with brake. Possibility to insert over the WC with standard dimensions. Demountable backrest for space saving during the transport.