



## 12.15XXL COMODA



**Categoria Prodotto:** [Comode](#), [Extra Large Solutions](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/12-15xxl-comoda/>

### Descrizione Prodotto

Aid with reinforced structure destined to very heavy persons.

Rigid reinforced frame made of atoxic powder varnished steel shapes. Padded seat

covered in leather like material very comfortable. The depth of the seat can be adapted to persons tall from cm 160 to 195. Available with or w/o WC. Backrest reclining until 60° and the articulation is reciprocal with the movement of inclination of the body.

Also the legrests are reclining, adjustable in flexion and extension and articulate reciprocally with the knee tibials. Removable armrests for easy side transfer and to approach tables. The rear wheels are Ø 600, pneumatic, reinforced with handrim and parking brakes.

**The front wheels are Ø 200x50 pneumatic assembled on ball bearings.**  
**Available seat sizes: cm 65-70-75. Max load 220 kg.**