



12.17N SERVO-ASSISTED JUNIOR VARIO PLUS



Categoria Prodotto: [Junior Line](#)

URL Prodotto:

<https://vassilli.it/en/prodotto/12-17n-servo-assisted-junior-vario-plus/>

Descrizione Prodotto

Suitable for childrens from 8 to 15 years

Indicated for users of limited mobility who need of continue differentiations of the position to compensate difficulties respiratory, post surgery, circulatory, articular and/or of physical resistance imposed by distonies and/or spasms. The frame is rigid made of atoxic powder varnished steel shapes. Seat adjustable in depth with reciprocal articulations. Armrests removable and adjustable in height. Footrests adjustable and removable. Tilting until 45° by gas pump, backrest inclination until 60° by gas pump. Padded contoured seat fixed with Velcro on tensioning bottom. Automatic wrapping backrest. Padded semiinvasive hadrest adjustable in height and depth on rotating base. **Available seat width cm 36. Max load kg 75.**