



12.17N XL VARIO PLUS



Categoria Prodotto: [Extra Large Solutions](#), [Comode](#)

URL Prodotto: <https://vassilli.it/en/prodotto/12-17n-xl-vario-plus/>

Descrizione Prodotto

Suitable for users with limited mobility which need continuous changements of the position to compensate the breathy difficulties, post surgery, circulatory, articular

and/or of physical strenght imposed by dystonia and/or spasms.

Rigid frame suitably reinforced, made of atoxic powder varnished steel shapes.

Seating adjustable in depth by reciprocal articulations. Armrests removable and adjustable in height. Footplates adjustable and removable. Tilting until 45° by gas pump, backrest inclination until 60° by gas pump. Padded sector seat fixed on the bottom with fabric removable. Contoured wrapping backrest on velcro. Padded headrest semi-invasive adjustable in height and depth on rotating base.

Complete with antitipping wheels. Available seat sizes: cm 50-55-60. Max load 160 kg.