



12.70N 4LIFE Junior



Categoria Prodotto: [Junior Line](#)

URL Prodotto: <https://vassilli.it/en/prodotto/12-70n-4life-junior/>

Descrizione Prodotto

Suitable for childrens from 8 to 15 years

Aluminium frame compact and reducible for the transport. Padded backrest, seating and backrest bottom tensioning adjustable. Backrest adjustable in inclination until 12° and in height until + 8 cm, completely forwards folding for space saving during the transport. Seat forwards adjustable in depth until + 6 cm

and in inclination until 12°. Side panels made of composite materials, adjustable in height and removable. Footrests made of duraluminium, adjustable in height and angle, with padded calfrest.

Front puncture proof wheels Ø 100 mm on duraluminium forks with holes for the adjustment of Ø front wheels; rear light alloy wheels Ø 600 mm, pneumatic and with quick release. Side structure for support of the rear wheels with locations for the changement of the depth; system for the adjustment of the convergence and of the camber from 0° to 8° of the rear wheels on request.

Seat adjustment in height from 44 to 55 cm with tilt +14° to -4°. **Available size cm 36 junior. Max Load kg 100. Registered model.**

With this product our Company has decided to engage itself to satisfy all the requests of users active on daily life who wants to have a product compact, adjustable, light and high functional.