



## 17.57N EVOLUTION ACTIVA



**Categoria Prodotto:** [Mobility](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/17-57n-evolution-activa/>

### Descrizione Prodotto

**Folding wheelchair with frame made of aluminium** and composite materials. Total weight kg 14. Separate removable footrests for space saving and for wheelchair folding. Fitted up with big wheels Ø 600 pneus with quick release and

front wheels Ø 120. Backrest with pushing handgrips adjustable in height.

Possibility to adjust the seating depth of the wheelchair until + 8 cm. Moving of the wheel axles and of the perpendicularly placing of the vertical axle of the front wheels. Camber of 5° on the rear wheels.

**Loading until 120 kg. Available with seat width 38, 40, 42, 44, 46, 48 cm.**

**Available sizes: 38 - 40 - 42 - 44 - 46 - 48 cm. Max load kg 120. Registered model.**