



## 17.57N JUNIOR EVOLUTION ACTIVA



**Categoria Prodotto:** [Junior Line](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/17-57n-junior-evolution-activa/>

### Descrizione Prodotto

Suitable for childrens from 8 to 15 years

Folding wheelchair with frame made of aluminium and composite materials. Total

weight kg 14. Separate removable footrests for space saving and for wheelchair folding.

Fitted up with big wheels Ø 560 pneus with quick release and front wheels Ø 120.

Backrest with pushing handgrips adjustable in height. Possibility to adjust the seating depth of the wheelchair until + 8 cm.

Moving of the wheel axles and of the perpendicularly placing of the vertical axle of the front wheels. Camber of 2° on the rear wheels.

**Loading until 100 kg. Available with seat width 36 cm. Registered model.**