



## 17.60XXL EVOLUTION



**Categoria Prodotto:** [Extra Large Solutions](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/17-60xxl-evolution/>

### Descrizione Prodotto

Folding reinforced frame with double cross made of steel varnished with powders, adjustment system allowing to increase the seat depth until 8 cm.

Synthetic loading fabrics padded, atoxic and hypoallergenic, washable and easy to be removed for replacement or sanitization. Backrest adjustable in height + 8 cm. Upsetting armrests with front lockage equipped with side containing protections. Rear pneumatic reinforced wheels Ø 600 mm, assembled on ball bearings and equipped with parking brake.

Support of rear wheels setted up for the moving of the pins. Separate adjustable legrests with one piece footplate. Revolving and removable footplate supports with quick release unlockage, complete with removable legrest strap. Fork support equipped with device for the adjustment of the vertical rotation axle.

Front directional pneumatic wheels Ø 200x50 mm. Rear handles with handgrips for the attendant. Rear device for easy step climbing.

**Available seating width: cm 65-70-75. Max Load kg 220.**