



17.68N EVOLUTION LIGHT



Categoria Prodotto: [Mobility](#)

URL Prodotto: <https://vassilli.it/en/prodotto/17-68n-evolution-light/>

Descrizione Prodotto

17.68N EVOLUTION LIGHT

Folding wheelchair with aluminium frame. Total weight 15 kg. Varnished structure. Footrests detachable, revolving and adjustable in height. Desk upsetting padded armrests. Tensioning backrest adjustable in height. Fitted up with rear wheels Ø

600 pneus with quick release and front wheels Ø 180. Possibility to change the load by shifting the wheels axle centre.

Available with seat width 40, 42, 44, 46, 48, 50, 52, 55 cm. Loading until 120 kg. Registered model