



## 17.68N JUNIOR EVOLUTION LIGHT



**Categoria Prodotto:** [Junior Line](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/17-68n-junior-evolution-light/>

### Descrizione Prodotto

Suitable for childrens from 8 to 15 years

Folding wheelchair with aluminium frame. Total weight 15 kg. Varnished structure.

Footrests detachable, revolving and adjustable in height. Desk upsetting padded armrests.

Backrest adjustable in height. Fitted up with rear wheels Ø 560 pneus with quick release and front wheels Ø 150. Possibility to change the load by shifting the wheels axle centre.

**Loading until 100 kg. Available with seat width 36 cm Junior. Registered model.**