



17.70N EVOLUTION ACTIVA COMPACT JUNIOR



Categoria Prodotto: [Junior Line](#)

URL Prodotto:

<https://vassilli.it/en/prodotto/17-70n-evolution-activa-compact-junior/>

Descrizione Prodotto

Suitable for childrens from 8 to 15 years

Evolution Activa Compact is a functional extension of the basic elements of our Evolution line. With the new Compact line, we achieved the maximum reduction of the overall dimensions, while keeping an adjustable seat which can suit also users of medium/large size.

Folding frame with double cross of duraluminium and composite materials.

Possibility of changing the seat depth. Possibility of changing the arrangement and the distribution of the loads by means of the slots on the big and small wheels.

Equipped with height adjustment. Onepiece footrest opening and adjustable in height and inclination of the foot. Padded seat and backrest with tensioning bottoms. Clothes protection side panels. Possibility of changing the height of the backrest. Parking brake. Adjustable handles for attendant.

Detachable rear wheels Ø 560 pneus and front ones puncture proof Ø 120. Bags for belongings on the backrest. Loading until 100kg.

Available with seat width 36 cm. Camber of 2° on the rear wheels.

Registered model.