



## 17.82N EVOLUTION ELECTRONIC JUNIOR



**Categoria Prodotto:** [Junior Line](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/17-82n-evolution-electronic-junior/>

### Descrizione Prodotto

Suitable for childrens from 8 to 15 years

Electric wheelchair with folding frame, of aluminium. Backrest and handgrips adjustable in height. Footrests adjustable, detachable, opening and retractable. Desk upsetting armrests.

Rear wheels Ø 300 mm and front wheels Ø 200 mm. Two motors of 160 Watts with gearbox, unlocking device (free wheels). Two sealed batteries of 42Ah maintenance free, only charging required and detachable. Programmable drive. Controller fitted up with joystick. Electromagnetic brakes.

The seating is provided with a system to increase the depth up to 8 cm. Loading until 120 kg. Autonomy of 35 km. Maximum speed programmable until 7,5 km/h.

**Available with seat ca 36 junior. Registered model.**