



## 17.82N EVOLUTION ELECTRONICS



**Categoria Prodotto:** [Mobility](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/17-82n-evolution-electronics/>

### Descrizione Prodotto

**Electric wheelchair with folding** varnished steel frame. Backrest and handgrips adjustable in height. Footrests adjustable, detachable, opening and retractable.

Desk upsetting armrests.

Rear wheels Ø 300 mm and front wheels Ø 200 mm.

Two motors of 180 Watts with gearbox, unlocking device (free wheels). Two sealed batteries of 4Ah maintenance free, only charging required and detachable.

**Programmable drive.** Controller fitted up with joystick. Electromagnetic brakes.

The seating is provided with a system to increase the depth up to 8 cm. Loading until 120 kg. **Autonomy of 35 km.** Maximum speed programmable until 7,5 km/h.

**Available with seat 42, 46 cm. Registered model.**