



18.64L HI-LO JUNIOR MANUAL GAS



Categoria Prodotto: [Junior Line](#)

URL Prodotto: <https://vassilli.it/en/prodotto/18-64l-hi-lo-junior-manual-gas/>

Descrizione Prodotto

Suitable for childrens from 8 to 15 years

To take the vertical position is not only a widening of autonomy, but it is also a therapeutical benefit. The verticalization allows to prevent many diseases connected to the motionlessness like the raise of decubitus, bad circulation, various pathologies of the muscular skeletal apparatus. The vertical movement is reciprocal to the articulations of the user
Verticalizing with levers and gas.

The new HI.LO is a technical aid very active easy to be used. The drawing is harmonious and easy to be customized thanks to the adjustments of the seating and backrest and setting of the wheels. It allows a distribution of the load reducing the effort for the moving.

The seating and the backrest are equipped with tensioning bottom covered with removable cushions, the depth can be adjusted; also the backrest in adjustable in height and upsetting forwards. Adjustable and removable footrests, side panels, chest padded strap adjustable and with quick release opening. Side levers which acts for opening the gas pump. Knee restrain adjustable and removable.

Warning for the supply of this wheelchair is necessary to have the dimensions of the user from seated position and also

his weight. **6 Available in seating width: cm 36 junior.**