



## 18.64MB HI-LO BAMBY



**Categoria Prodotto:** [Junior Line](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/18-64mb-hi-lo-bamby/>

### Descrizione Prodotto

Suitable for childrens from 3 to 6 years and from 6 to 10 years

Frame made of duraluminium varnished with epoxy powders. Wheelchair with

system of verticalization by gas pump with lever control on the back of the backrest operating by the attendant. Padded seat and backrest covered in washable material and easily removable for the sanitization, seating and the backrest adjustable in tension.

Backrest adjustable in height and inclination. Adjustment of the seating depth of + 10 cm. Padded chest strap which can be opened by quick release. Side panels. Push handles for the attendant adjustable in height. Removable footrest with independent aluminium footplates, adjustable in height, depth and inclination. Removable padded knee restrain applied to the footrest adjustable in height width, depth and inclination. Front directional puncture proof wheels Ø 100 mm. Rear pneumatic wheels Ø 560 mm with handrim for the selfpropelling removable and with quick release. Footrest support and knee restrain easily removable.

**Available seating width: from cm 28 to 34. Max Load kg 75.**