



18.64MBP BAMBY PLUS



Categoria Prodotto: [Junior Line](#)

URL Prodotto: <https://vassilli.it/en/prodotto/18-64mbp-bamby-plus/>

Descrizione Prodotto

Suitable for childrens from 8 to 15 years

Compact duraluminum frame varnished with epoxy powders with seating

adjustable in depth and height. Verticalizing seating adjustable in depth and width, tilting until 32° by independent gas pumps with mechanical lever controls on the handgrips for the attendant.

Backrest adjustable in height and inclination by gas pump with handles for the attendant adjustable in height. Side panels for pelvis containing. Removable footrest with independent aluminium footplates, adjustable in height, depth and inclination.

Removable padded knee restrain applied to the footrest adjustable in height, width, depth and inclination. Footrest support and knee restrain easily removable. Padded PU armrests with front lockage adjustable in height, width and depth. Padded half containing headrest, fully adjustable in depth, height and inclination.

Padded chest strap, applied to the backrest, adjustable in width and opening by quick release. Directional front puncture proof wheels Ø 100 mm.

Rear pneumatic wheels Ø 320 mm with possibility of placing on various positions to give a customized arrangement of seating and movement. Pressure brake on the rear wheels.

Available seating width: cm 28 - 34. Max Load kg 75. Seat and backrest paddings NOT included.