



18.70 HI-LO VARIO Junior



Categoria Prodotto: [Junior Line](#)

URL Prodotto: <https://vassilli.it/en/prodotto/18-70-hi-lo-vario-junior/>

Descrizione Prodotto

Suitable for childrens from 8 to 15 years

The new electronic wheelchair HI-LO VARIO satisfy the priority need in the assistance and care of patients having deep neurologic and motor deseases

because it improves the circulation and the functioning of the gastrointestinal apparatus, reduces the arising of osteoporosis, improve the trofism muscle-skeletal, reduces the risk of arising of ulcers and/or decubitus sores.

The daily use of an aid such HI-LO VARIO carry out a “preventive” and “rehabilitative” function due to the technical features highly specific and modular. The seat and the backrest are supplied with tensioning belts that allows the adaptation to the structure of the user, and covered by a padded cushion covered by an elastic tissue, the articulations follow reciprocally the movements of the user, the depth and the width of the seating are adjustable to be adapted to the different sizes of the users; the front wheels opens automatically in order to increase the stability.

All the movements are activated by a joystick by which is possible to: - adjust the seating height; - tilt the seating; - verticalize the seating; - recline the backrest and the footrests. The knee restrains are adjustable and removable by quick release device, the headrest is adjustable in height and depth, the armrests are upsettable and adjustable in angle with safety front closing.

Autonomy until 40 km on flat surfaces. Available in size: 36 cm Junior