



18.80 JUNIOR NEW 4x4



Categoria Prodotto: [Junior Line](#)

URL Prodotto: <https://vassilli.it/en/prodotto/18-80-junior-new-4x4/>

Descrizione Prodotto

-- OUT OF PRODUCTION --

Suitable for childrens from 8 to 15 years

It is a product designed by the Centre of Research Tecnothon of Foundation Telethon in cooperation with Vassilli company. 4x4 is a wheelchair very easy to be used, which adapts to every kind of ground, also sandy due to the integral drive, by its reduced dimensions it suits very well also on home areas.

It is a technical aid having high functional features, the design is charming and modern, the structure is made of duralluminium and steel shapes, the front wheels drive support is articulated; the seating is adjustable in height from cm 43 to 55 by electric actuator, adjustment in depth until + 12 cm.

The backrest is electrically reclining with functions directly on joystick, armrests adjustable and removable, adjustable headrest, footrests adjustable and removable.

Available with seating width of cm 36 junior.