



18.82N ELECTRONIC



Categoria Prodotto: [Mobility](#)

URL Prodotto: <https://vassilli.it/en/prodotto/18-82n-electronic/>

Descrizione Prodotto

POWERWHEELCHAIR

The new 18.82N is a very versatile and functional product. With the technical solutions with which is equipped the wheelchair it is possible to adjust both the seating depth from cm 40 to 52, the width of + 7 cm on all the basic sizes. The

armrests are adjustable in height and upsetting.

The seating and the backrest are composed of tensioning belts and a padded detachable cushion. The backrest can be longered with angle adjustable until 12° and is foldable forwards for space saving for transport. The footrests are revolving, adjustable and removable.

It is possible to adjust the seating angle on the telescopic column under the seat until 16°. The wheels are pneumatic with increased section. The electronic components have a power until 50 Ah, the batteries are 42 Ah maintenance free; it is equipped with electronic battery charger. Autonomy until 40 Km in standard conditions.

User load capacity until 120 Kg. Available with seating of 42 and 46 cm.