



18.82N ELETTRONICA JUNIOR



Categoria Prodotto: [Junior Line](#)

URL Prodotto: <https://vassilli.it/en/prodotto/18-82n-elettronica-junior/>

Descrizione Prodotto

Suitable for childrens from 8 to 15 years

The new 18.82N is a very versatile and functional product. With the technical solutions with which is equipped the wheelchair it is possible to adjust both the seating depth from cm 40 to 52, the width of + 7 cm on all the basic sizes.

The armrests are adjustable in height and upsetting. The seating and the backrest

are composed of tensioning belts and a padded detachable cushion. The backrest can be longered with angle adjustable until 20° and is foldable forwards for space saving for transport. The footrests are revolving, adjustable and removable. It is possible to adjust the seating angle on the telescopic column under the seat until 16°.

The wheels are pneumatic with increased section. The electronic components have a power until 50 Ah, the batteries are 42 Ah maintenance free; it is equipped with electronic battery charger. Autonomy until 40 Km in standard conditions.

User load capacity until 120 Kg. Size of the seating cm 36