



## 19.98N JUNIOR NEW SPACE 1



**Categoria Prodotto:** [Junior Line](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/19-98n-junior-new-space-1/>

### Descrizione Prodotto

Suitable for childrens from 8 to 15 years

Electronic wheelchair mainly for outdoor use which allows to move easily also indoor thank to the compact dimensions. The design is modern and charming, it includes technical and innovative solutions.

The seat and the backrest are adjustable in depth and height by tensioning bottoms and padded with removable cushion lycra covered. Padded armrests adjustable in height, removable. Removable side panels adjustable in width until 7 cm on each basic frame. The backrest can be upsetted forwards and has an adjustable angle until 12°, the angle of the seat can be adjusted until 10° backwards. Adjustable and removable footrests.

Pneumatic increased wheels, 2 geared motors of 350 Watts, 2 batteries 65Ah free maintenance. Joystick and electronic controller of 80Ah. Set of lights and indicators. Equipped with 8Ah battery charger. 5 levels of speed until max 12 km. Autonomy in standard conditions until 5 hours. Load capacity of user until 120 kg.

**Available in seat width cm 36 junior.**