



## 19.98XXL 280 kg NEW SPACE



**Categoria Prodotto:** [Extra Large Solutions](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/19-98xxl-280-kg-new-space/>

### Descrizione Prodotto

This aid has been made to support bariatric users weighting until 280 kg. The structure is made of reinforced steel varnished with epoxy powders. Motors are

enhanced 500W and there are 2 free maintenance batteries of 100Ah each one, the electronic system is P&G of 120Ah.

The seat frame is adjustable in width and depth to adapt to an eventual increase of the user's dimensions. Pneumatic increased wheels, footrests removable and adjustable in height, seat and backrest on rigid bottom with PU high density foam cushion 6 cm of height covered in breathable fabric and/or in waterrepellent PU, the seat depth can be adjusted from cm 45 to 60, the backrest is predisposed for angle adjustment from 0° to 15° which can be set up on fix points each 3°, the armrests are adjustable in height with side protections, they can be upsetting for the side transfer and the approach to tables and desks.

**The aid is available in three basic seat width: cm 65 - 70 - 75 - max load kg 280.**