



20.12 Siliconed fibre mattress



Categoria Prodotto: [Beds and Homecare Aids](#)

URL Prodotto: <https://vassilli.it/en/prodotto/20-12-siliconed-fibre-mattress/>

Descrizione Prodotto

Padding made of siliconed fibre contained in easily detachable cylinders.

The fibre ensure to the patient a good distribution of the resting pressures. This mattress was expressly designed to allow also the complete pressure unloading detaching one or more cylinders from the most critical areas.

The external cover, with zip throughout one side is separated in 4 sections where the 13 cylinders are located.

Dim. cm 190x83x12h. Requirements following Rule ISO EN 12182.