




20.14 Composed mattress



Categoria Prodotto: [Beds and Homecare Aids](#)

URL Prodotto: <https://vassilli.it/en/prodotto/20-14-composed-mattress/>

Descrizione Prodotto

Lower and external part made of high density cold foam.

Partition of the mattress in 2 parts of 5 cm of height. Inside there is a surface composed of detachable squares made of viscoelastic material for the differentiated discharged of the pressure and a surface with squares with exhaust below to allow the inner ventilation.

Possibility to insert the siliconed fibre section mattress.

Dim. cm 195x90x15h.

Requirements following Rule ISO EN 12182.