



20.20 Air mattress with flating elements



Categoria Prodotto: [Beds and Homecare Aids](#)

URL Prodotto: <https://vassilli.it/en/prodotto/20-20-air-mattress-with-flating-elements/>

Descrizione Prodotto

Possibility to adjust the pressure following the user weight.

The elements of the mattress of nylon + PU are fixed by buttons easily interchangeable in case of breakage and washable. The special pump particularly silent allows an optimal adjustment of the pressure by the handle on the front of the device.

Two hooks allow to hang the pump to the bed edge.

Dim. cm 200X90X11,5h. Requirements following Rule EN 12182 - CEI UNI EN 60601-2-52.