



20.22 Air bubble mattress



Categoria Prodotto: [Beds and Homecare Aids](#)

URL Prodotto: <https://vassilli.it/en/prodotto/20-22-air-bubble-mattress/>

Descrizione Prodotto

The alternate of the points of contact with the skin effects a soft massage which allows a good blood flow and a good oxygenation of the tissues, preventing the arise of decubitus ulcers.

Easy installation, use and maintenance. The mattress is made of Medical PVC, particularly comfortable, it is fixed to the mattress by flaps. All the components of the mattress are washable.

Dim. cm 200X90X6h.

Requirements following Rule ISO EN 12182 - CEI UNI EN 60601-2-52