



## 400.20 WC stabilizing uprised seat



**Categoria Prodotto:** [IHS](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/400-20-wc-stabilizing-uprised-seat/>

### Descrizione Prodotto

Structure made of steel varnished tubular with vinyl seat and plastic bowl with antislip rubber tips, adjustable in height. It includes various functions together: uprise for the WC, support for the WC, commode chair. All the material used for the manufacture are hypoallergenic and atoxic.

Load capacity kg. 100.

Aid in conformity to the Rule EN 12182.