



## 400.30 Bathtub revolving seat



**Categoria Prodotto:** [IHS](#)

**URL Prodotto:** <https://vassilli.it/en/prodotto/400-30-bathtub-revolving-seat/>

### Descrizione Prodotto

Structure made of tubular inox steel with comfortable vinyl seat and rubber antislip tips for bathtub support. The seat and the armrests rotate at 360° in steps of 90° and their dimensions are compact inside the external dimensions of the frame. In the lower part of the frame, on both sides are placed 2 stabilizers screw adjustable (width max 67 cm - width min. 51 cm). All the materials used for the manufacture are hypoallergenic and atoxic.

Load capacity kg. 100.

Aid in conformity with the Rule EN 1218.